



Product Code	27300
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279273007

Description

GIMA Mechanical baby scale

High quality sturdy mechanical baby scale for accurate and sensitive weighing.
Supplied with curved weighing tray, aluminium levers and pins, stable casing for safe weighing.
User manual in GB, FR, IT, ES, PT, DE, RO, GR.

Technical Specifications

Capacity: 20 kg
Sensitivity: 10 g
Dimensions: 570x330x285 mm
Weight: 6.5 kg