

# SOEHNLE 8352 MULTINA DIGITAL BABY AND CHILD SCALE



Product Code	27306
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279273069
---------	---------------

PARAF	903170468
-------	-----------

## Description

### SOEHNLE 8352 MULTINA DIGITAL BABY AND CHILD SCALE

Well-designed, easy to handle plastic baby scale.

- suitable for children up to 20 kg removing baby tray
- easy to read LCD-display with 13 mm figures
- simple operation with one button
- tare function to calibrate for baby blankets/textiles
- fully-electronic scale: no wear and tear parts, no maintenance
- hold function: weight value will be held for 10 seconds after unloading
- automatic power-off

Made in Germany.

Manual in: DE, GB, FR, IT, ES, NL, DK, SE, HU, FI, PT.

## Technical Specifications

Capacity: 20 kg  
Sensitivity: 10 g  
Dimensions: 601x346x92 mm  
Weight: 3.5 kg  
Battery: 9 V  
Warranty: 2 years