



Product Code	27312
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279273120
PARAF	903605184

Description

BABY AND CHILD DIGITAL SCALE

Multifunctional scale for babies and children with removable light plastic tray.
 It can be used with tray for babies or without for children.
 Weight difference display function with memory: you can see the gain of weight of the baby after feeding.
 Precise scale with easy-to-read LCD display.
 Auto zero, auto tare, auto hold function.
 Automatic switch off and low battery indicator.
 Switchable to kg/lb/st.

Multilanguage manual and box in: GB, FR, IT, ES, PT, PL, DE, RO, CZ, HR, HU, SI.

Technical Specifications

Capacity: 20 kg
 Sensitivity: 5 g
 Tray dimensions: 525 x 305 x 63-95 mm
 Weight: 2 kg
 Pack dimension: 540 x 300 x 110 mm
 Battery operated: 2 x CR2032 lithium batteries