



| | |
|---------------|-------------------|
| Product Code | 27313 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

| | |
|---------|---------------|
| UK-REP | Not required |
| CH-REP | Not required |
| EAN/UPC | 8023279273137 |
| PARAF | 903118053 |

Description

GIMA ELECTRONIC BABY SCALE

Handy, reliable and well-designed electronic baby scale.

Provided with large LCD display 66x28 mm and kg/lb switch.

Supplied with height tape meter and washable sheet.

Multilingual user manual and gift box: GB, FR, IT, ES, PT, DE, Arabic

LCD digit height: 18 mm

Technical Specifications

Capacity: 20 kg / 44 lbs
 Sensitivity: 10 g
 Dimensions: 595x340x90 mm
 Weight: 4 kg
 Battery operated: 4x1.5 V AA
 Auto-off function