



Product Code	27317
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No
RDM (NSIS)	1412922
CND	V0399
EAN/UPC	4012030008945

## Description

### SECA 213 MOBILE STADIOMETER

Mobile stadiometer for measuring height, simple and easy to transport with integrated handle.  
 The measuring rod can be dismantled into several pieces and set up easily and quickly.  
 The spacer keeps the rod straight and stable without any fittings.  
 Flexible position due to large floor plate which ensures stability.  
 It guarantees precise and clearly visible reading of results.  
 Suitable for adults and children and for measuring patients lying down.

Multilingual manual in: GB, FR, IT, ES, PT, DE, PL, NL, GR, DK, SE, NO, FI.

## Technical Specifications

Measuring range: 20-205 cm  
 Graduation: 1 mm  
 Size (open): 337x590xh 2,130 mm  
 Size (folded): 337x177x624 mm  
 Weight: 2.4 kg