

PLICOMETRO - SKINFOLD CALIPER 0-40 mm - mechanical



| | |
|---------------|----------------|
| Product Code | 27320 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I M |

| | |
|------------|---------------|
| UK-REP | No |
| CH-REP | No |
| RDM (NSIS) | 2520432 |
| CND | V030299 |
| EAN/UPC | 8001066702008 |
| PARAF | 977811987 |

Description

PLICOMETRO - skinfold caliper 0-40 mm - mechanical

This high precision instrument measures the thickness of a fold of skin with its underlying layer of fat giving physicians the situation of the subcutaneous tissues of the patient.

Possibility of re-calibration

Produced to meet International Standards, the Gima "Skinfold caliper" has springs which exert a pressure of 10 g/mm² on the skinfold and an accurate scale which measures thickness in mm. With possibility of re-calibration of the pressure it solves the major problem of traditional skinfold caliper that is reduction of the tension of the spring in the time.

Easy to use and read

With double reading scale, it is supplied with picture showing advised positions for measurements. Tables will help obtain correct % of body fat from measured data and make a correct evaluation of the patient.

Stylish and user-friendly

A special nylon finish makes the device feel skin, avoiding trouble of contact with metal. Made of metal and sold in an attractive wooden box (19x36xh 6 cm).

User manual in GB, FR, IT, ES, DE.

Made in Italy.