



| | |
|---------------|----------------|
| Product Code | 27326 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I M |

| | |
|--------|----|
| UK-REP | No |
| CH-REP | No |

| | |
|------------|---------------|
| RDM (NSIS) | 2520420 |
| CND | V030201 |
| EAN/UPC | 8001066104550 |

Description

BABY MEASURE "CALIBRO"

Made of oxidized aluminium, metric graduation 0-900 mm, division 5 mm.

Certified by Italian "weights and measure institute".

Weight: 300 g.

Multilingual manual: GB, FR, IT, ES, PT, DE, PL, SE, RO, GR.

Made in Italy.