



Product Code	27340
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279273403
PARAF	903174946

## Description

GONIOMETER - 205x45 mm

With circular scale used to check the size of hip joints and the angle of the femoral neck on X-rays.  
Personalisation with your logo for minimum 1,000 pcs.

Used by orthopaedic surgeons to measure the range of movement of joints: elbows, knees and fingers.  
Approximate measurement of values.

Made in the UK.