

# FAT 1 SKINFOLD CALIPER



Product Code	27344
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	Yes

RDM (NSIS)	787496
CND	V030299
EAN/UPC	8023279273441
PARAF	922289677
GMDN	32262

## Description

### FAT-1 SKINFOLD CALIPER

Economic but precise and accurate skinfold caliper for professional use.

Measures the thickness of skin folds over the body, based on the characteristics developed by Tanner and Whitehouse:

- pressure on the fold: approx 10 g/mm<sup>2</sup>
- contact surface on the fold: approx 6x11 mm

It allows a correct evaluation of the nutritional status and the sectoral distribution of adipose tissue.

Raise a fold of skin with your fingers and measure its thickness by the forceps (max opening 4 cm): on the top edge a scale in millimeters and an arrow measure the skin fold (triceps, subscapular muscles, etc...) with an accuracy of about  $\pm 2$  mm.

The coloured bands placed near the scale allow a simultaneous comparison of the measured value with referential percentages related to men and women.

Advised in dietetics, cardiology, sports medicine, aesthetic medicine, gynaecology, pediatrics, ecc...

Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, LT, PL, CZ, SK, HU, RO, BG, GR, Arabic.

Made in Italy.