



|               |                   |
|---------------|-------------------|
| Product Code  | 27350             |
| Unit of sale  | 1 pc              |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|         |               |
|---------|---------------|
| EAN/UPC | 8023279273502 |
| PARAF   | 984952717     |

## Description

PROFESSIONAL BODY TAPE with pediatric gauge box

Professional body tape measure for adult and pediatric use.

Tape is coloured on both sides to help to identify metabolic risk in abdominal measure for male and female.

Length of tape: 8-150 cm.

Thanks to its rounded shape, it can be used also as a head measure tape.

Packed in single box.

Body tape with pediatric gauge box. Specific tape for pediatric use.

Supplied in a box with a ruler to compare the measured head circumference with the standard percentiles of children's growth from 0 to 36 months.

Suitable for both sexes.

Approximate measurement of values.

Multilingual manual: GB, FR, IT, DE, ES, GR, PT, Arabic.