



Product Code	27353
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	762096
CND	V9099
EAN/UPC	8023279273533
PARAF	960149211

## Description

### LED PODOSCOPE

Device for the plantar analysis of the feet. It has a lacquered wooden structure and a crystal top with an inclined mirror below it for a real view of the sole of the feet, which highlights the points of greatest and least load. Double lateral illumination with high level effect lamps. Max load 200 kg. Strong, reliable and professional product.

Working surface 32x36 cm  
 Size: 46x55xh 33 cm  
 Weight: 15.5 kg  
 Power: 100-240 V, 50/60 Hz

Made in Italy.