

"ARMAGNAC" ILLUMINATED CHART 5 m



Product Code	27380
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	102692
CND	Z12120199
EAN/UPC	8023279273809
PARAF	903120499
GMDN	36380

Description

Armagnac Optometric chart

- Luminous model
- Reading distance: 5 meter

GIMA luminous optotypes have been designed following clinical recommendations for an extra white light to give an accurate and precise reading of any detail.

Manufactured with acrylic screen, and according to safety norms (fuses, Shuko plug, green light On/Off,...).

Made in Italy.

Technical Specifications

- Luminous model
- Reading distance: 5 meter
- Viewing area: 24 x 62 cm
- External dimensions: 29.5 x 67.5 x 12 cm
- Weight: 6 Kg