



Product Code	27387
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	1433462
CND	V9099
EAN/UPC	8023279273878

Description

Device for analysing posture in the frontal, rear and lateral planes.
 The patient stands with feet in the positions indicated on the platform.
 It includes a platform, elastic grid, two aluminum side bars, measuring indicators with plumb line for postural reference (Barré's vertical evaluation) and adjustable mirror on top.
 The image reflected in this upper mirror allows observation of any rotation of the shoulders.

Manual: GB, IT. On request ES.

Supplied with a set of three horizontal and two vertical elastics that are respectively adjustable in height and width by sliding. These elastics can be used as reference lines in the observation of the patient's upright posture, and in the search for asymmetries and deviations.

Made in Italy.

Size: 60 x 66 x h 220 cm
 Capacity: 200 kg