

CLASSIC EXAMINATION COUCH - with hole



Product Code	27405
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2537687
CND	V080603
EAN/UPC	8023279274059
GTIN	08023279274059
GMDN	38458

Description

CLASSIC EXAMINATION COUCH - chromed - blue- backrest with hole (hole 23x12 cm)

Chromed mainframe and legs constructed in chrome-plated steel.

Legs fitted with plastic feet.

The two-section top has resilient foam padding covered with washable and fire-proof plastic material.

Backrest fully adjustable.

Size: 180x57xh 73 cm.

Maximum load 120 kg.

Delivered in kit form (198x17x64 cm).

Instructions: GB, FR, IT, ES, PT, DE, GR.

Fireproof, simulated leather, washable upholstery in two sections

Made in Italy.