

SUITCASE COUCH



Product Code	27407
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2537691
CND	V080603
EAN/UPC	8023279274073
PARAF	903112922
GMDN	38458

Description

SUITCASE COUCH - black

Chromed steel and folding legs.
Upholstered covering in washable material.
Resin padded.
Max load: 100 kg.
Weight: 16.5 kg
Supplied with a velvet bag for easy transport.
Packing: 95x18x73.5 cm
Manual: GB, IT, NL.

Couch size:
187.5x64xh 71 cm (open)
93x68xh 14 cm (folded)

Mattress thickness 3 cm

Made in Italy.