



Product Code 27416

Unit of sale 1 pc

Minimum order 1

Type No medical device

Class

UK-REP Not required

CH-REP Not required

EAN/UPC 8023279274165

## Description

WOODEN COUCH - beech

MAXIMUM COUCH LOAD: 200 KG

Comfortable couches, ideal for physical therapy, esthetics and wellness centres.

Superior quality wood.

Back rest adjustable in 13 positions.

Resilient foam padding covered with washable skay.

Padding in beige colour.

Size: 190x70xh 75 cm.

Delivered in kit form. Easy assembly, tool provided.

Made in Italy.