

# DELUXE TROLLEY with 2 drawers 29 x 40 X H 10.5 cm



Product Code	27422
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279274226
GTIN	08023279274226

## Description

DELUXE TROLLEY with 2 drawers

Size: 65x45xh 80 cm

Drawer size: ext. 58x40xh 10.5 cm, int. 51.5x38.5xh 8.5 cm.

Top quality trolleys made of lightweight anodized special alloy with detachable stainless steel shelves.

4 revolving rubber castors Ø 80 mm.

Max load (x shelf): 20 kg.

Instructions: GB, FR, IT, ES PT, DE, GR.

Delivered in kit form.

Made in Italy.