

MICROEMBOSSED GLUED 2 PLIES PEDIATRIC COUCH ROLL - 46m x 50cm



| | |
|---------------|-------------------|
| Product Code | 27423 |
| Unit of sale | 9 rolls |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

EAN/UPC 8017414001327

Description

MICROEMBOSSED GLUED 2 PLIES PEDIATRIC COUCH ROLL - 46m x 50cm - box of 9

Decorated medical rolls for pediatric use. Produced in pure cellulose, dermatologically tested.
Hypoallergenic inks.
Each box contains 9 rolls with 5 assorted themes.

Material: Microembossed glued
Colour: white+decoration
Ply: 2 ply
Weight: 2x17.5 g/m²
Length: 46 m
Height: 50 cm
Number of sheets: 135
Perforation: 34 cm
Rolls per box: 9