

MICROEMBOSSED GLUED 2 PLYS COUCH ROLL - 50m x 50cm



Product Code	27424
Unit of sale	9 rolls
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2816520
CND	T030502
EAN/UPC	8023279274240

Description

COUCH PAPER ROLLS

- Material: Microembossed glued
- Colour: white
- Ply: 2 ply
- Weight: 2 x 17.5 g/m²
- Length: 50 m
- Height: 50 cm
- Number of sheets: 132
- Perforation: 38 cm
- Rolls per box: 9

A complete range of hygienic and resistant couch paper manufactured in Italy in pure cellucotton or cellucotton/polythene. Each roll is individually packed in plastic bag for best protection. Atoxic. Available with standard 38 or 40 cm perforation or 190 cm perforation (full person size).