

GIMA 2 TROLLEY with guard-rail - small



Product Code 27429

Unit of sale 1 pc

Minimum order 1

Type No medical device

Class

UK-REP No

CH-REP No

EAN/UPC 8023279274295

PARAF 903605816

Description

GIMA 2 TROLLEY - small: 60x40xh 77 cm - weight 11 kg

Practical trolleys with aluminium body and handles.

Two white bilaminated shelves 18 mm thick, black PVC bordered, with guard-rail.

Two large lateral handles facilitate transport or lift up of the trolley.

4 revolving nylon castors Ø 80 mm.

Max load (x shelf): 40 kg.

Delivered in kit form.

Made in Italy.