


MICROEMBOSSED GLUED 2 PLYS COUCH ROLL - 80m x 59cm



Product Code	27432
Unit of sale	6 rolls
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279274325
---------	---------------

Description

COUCH PAPER ROLLS

- Material: Microembossed glued
- Colour: white
- Ply: 2 ply
- Weight: 2 x 17.5 g/m²
- Length: 80 m
- Height: 59 cm
- Number of sheets: 210
- Perforation: 38 cm
- Rolls per box: 6

A complete range of hygienic and resistant couch paper manufactured in Italy in pure cellucotton or cellucotton/polythene. Each roll is individually packed in plastic bag for best protection. Atoxic. Available with standard 38 or 40 cm perforation or 190 cm perforation (full person size).