



Product Code	27437
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279274370
PARAF	903605931
GTIN	08023279274370

Description

ANAESTHETICS TROLLEY

Made of satined AISI 304 stainless steel, with 6 drawers (internal 26x37xh 7 cm), a lockable vertical compartment and an open space. Top with 3.5 cm upstand. 4 revolving castors Ø 80 mm.

Size: 60x40xh 110 cm
Weight: 39 kg

Drawers are locked once vertical compartment is locked

Made in Italy.