



Product Code	27438
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279274387
PARAF	903113025

Description

TRIS TROLLEY

A practical trolley for ECG, Electrosurgical unit or any other device, with 4 nylon castors and 2 lateral handles for easy moving.

Size: 49x37xh 77 cm.

Shelf size: 40x30 cm.

Total load: 39 kg, max per shelf 13 kg

Colour: RAL 9002

Multilingual manual: GB, FR, IT, ES, PT, DE.

Delivered in kit form.

Revolving castors Ø 50 mm

Made in Italy.