

DELUXE TROLLEY



Product Code	27460
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8023279274608
PARAF	903113090

Description

DELUXE TROLLEY

Size: 65x45xh 80 cm
Top quality trolleys made of lightweight anodized special alloy with detachable stainless steel shelves.
4 revolving rubber castors Ø 80 mm.
Max load (x shelf): 20 kg.

Instructions: GB, FR, IT, ES PT, DE, GR.
Delivered in kit form.

Made in Italy.