



| | |
|---------------|-------------------|
| Product Code | 27460 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

| | |
|--------|----|
| UK-REP | No |
| CH-REP | No |

| | |
|---------|---------------|
| EAN/UPC | 8023279274608 |
| PARAF | 903113090 |

Description

DELUXE TROLLEY

Size: 65x45xh 80 cm

Top quality trolleys made of lightweight anodized special alloy with detachable stainless steel shelves.

4 revolving rubber castors Ø 80 mm.

Max load (x shelf): 20 kg.

Instructions: GB, FR, IT, ES PT, DE, GR.

Delivered in kit form.

Made in Italy.