

# OVERBED TABLE - moulded tray



Product Code	27469
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279274691
---------	---------------

## Description

OVERBED TABLE - moulded tray

Frame made of chrome plated steel tube.

Serving tray made of thermoformed anti-shock plastic material with retaining edges and pre-formed bottle, glass and cutlery holder.

Tray dimensions: 58 x 39 x h 5 cm

Adjustable height: 75 to 114 cm

Base with 4 Ø 50 mm castors.

Delivered in kit form. Carton size: 59 x 22 x 63 cm.

Made in Italy.