



Product Code	27481
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279274813
PARAF	983761824

## Description

### DRESSING TROLLEY - Standard

Colour: yellow

Size of trolley (without accessories): 82 x 64 x h 100 cm

Number of drawers: 4

Size of drawers:

- 60 x 48 x h 13 cm (3)

- 60 x 48 x h 27 cm (1)

These multifunctional professional trolleys have painted steel structure and external plastic finishing.

Provided with 4, 8 or 9 drawers, one lateral waste container (manual opening, capacity 15 l) and one extra lateral table.

Four swivelling rubber castors Ø 125 mm (2 with brake).

Equipped with instrument tray and upper dispenser (10 tilting bins).

Made in Italy.