

# RECORD HOLDERS TROLLEY - sliding top



Product Code	27490
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279274905
PARAF	983761836

## Description

### RECORD HOLDERS TROLLEY - sliding top

Record trolleys are designed to hold clinical reports in A3/A4 format or specific format.  
Frame in epoxy painted steel, working surface and panels in technopolimer plastic.  
Equipped with 4 Ø 125 mm swivelling rubber castors (2 with brake).  
Bumpers around the perimeter.

Size: 820x640xh 1,000 mm

Keeps 40 file racks (not included).

Made in Italy.