



| | |
|---------------|-------------------|
| Product Code | 27491 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

| | |
|---------|---------------|
| UK-REP | Not required |
| CH-REP | Not required |
| EAN/UPC | 8023279274912 |

Description

RECORD HOLDERS TROLLEY

Record trolleys are designed to hold clinical reports in A3/A4 format or specific format.
 Frame in epoxy painted steel, working surface and panels in technopolimer plastic.
 Equipped with 4 Ø 125 mm swivelling rubber castors (2 with brake).
 Bumpers around the perimeter.

Size: 820x640xh 1,000 mm

Keeps 25 file racks (not included).

Made in Italy.