

RECORD HOLDERS TROLLEY 2 large drawers



Product Code	27492
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279274929

Description

RECORD HOLDERS TROLLEY - 2 large drawers

Record trolleys are designed to hold clinical reports in A3/A4 format or specific format.
Frame in epoxy painted steel, working surface and panels in technopolimer plastic.
Equipped with 4 Ø 125 mm swivelling rubber castors (2 with brake).
Bumpers around the perimeter.

Size: 820x640xh 1,150 mm

Keeps 25 + 25 file racks (not included).

Made in Italy.