



Product Code	27501
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2537726
CND	V080603
EAN/UPC	8023279275018
PARAF	903141188
GTIN	08023279275018
GMDN	36065

Description

GYNAECOLOGY BED - standard leg holders

Strong chrome-plated structure. Rectangular steel tube. Adjustable head and leg sections. Simulated leather, black colour, washable upholstery, self extinguishing. Load: 110 kg. Complete with s/s bowl Ø 31 cm and leg-holders.

Size: 180x55xh 80 cm

Seat size: 55x55 cm

Delivered in kit form

Multilingual instructions: GB, FR, IT, ES, PT, DE, GR

Made in Italy.