

# 3-SECTION GYNAECOLOGICAL BED - blue



Product Code	27505
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2254549
CND	V0880
EAN/UPC	8055136920038

## Description

### 3-SECTION GYNAECOLOGICAL BED - blue

Strong and stylish structure in oval tube 50x25 mm, painted with epoxy powder, metallic grey colour. The three sections top has high density padding, covered with washable, fire proof, leatherette material Class 1 IM. Back-rest and leg-rest section are fully adjustable by self mechanical sytem. Gynaecological bed is supplied with 2 foam padded leg holders and stainless steel bowl Ø 32 cm. Maximum load: 150 kg - Weight: 36 kg.

Seat: 50x60 cm  
Size: 181x60xh 83 cm  
Delivered in kit form.

Made in Italy.