

HEIGHT ADJUSTABLE GYNAECOLOGICAL BED - blue



Product Code	27507
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2262252
CND	V08060199
EAN/UPC	8023279275070

Description

HEIGHT ADJUSTABLE GYN BED - blue

Strong and solid structure in oval tube 50x25x2 mm painted with epoxy powder, white colour. The three sections top has high density padding, covered with washable, fire proof, leatherette material Class 1 IM.

Fully adjustable

Height, up to 95 cm, is adjustable by electric engine (6000 N power) with foot control, head section is adjusted by a gas cylinder (400 N power), central section by a gas cylinder (600 N power) up to 65°.

This gyn bed is supplied with 2 quality foam leg holders, stainless steel tray and couch roll holder.

Maximum load: 150 kg.

Weight: 70 kg.

Stainless steel tray 370 x 270 x 50 mm

Size: 195 x 60 x h 48/95 cm

Made in Italy.