

DELUXE DELIVERY BED - with Trendelenburg



Product Code	27510
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2537726
CND	V080603
EAN/UPC	8023279275100
PARAF	903141190
GTIN	08023279275100
GMDN	36065

Description

DELUXE DELIVERY BED - quality leg holders
3 sections with independent controls and Trendelenburg position.

Strong chrome-plated structure. Rectangular steel tube. Adjustable head and leg sections. Simulated leather, black colour, washable upholstery, self extinguishing. Load: 110 kg. Complete with s/s bowl Ø 31 cm and leg-holders.

Size: 180x55xh 80 cm
Delivered in kit form

Multilingual instructions: GB, FR, IT, ES, PT, DE, GR

Made in Italy.