

# 3-SECTION DELIVERY BED WITH TRENDELENBURG



Product Code	27538
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	1229129
CND	V0899
EAN/UPC	8023279275384

## Description

### 3-SECTION DELIVERY BED - with Trendelenburg

- structure made of epoxy painted oval steel tube
- three sections with independent control and Trendelenburg position
- padding made of high density polyurethane foam, fireproof Class 1 IM
- covering made of washable, plastic material (leatherette) fireproof Class 1 IM
- adjustable rubber feet

Size: 181x60xh 80 cm

Carton size: 75x20xh 196 cm.

Maximum load: 150 kg.

Net weight: 43 kg.

Supplied with a pair of height-adjustable, padded, anatomic leg-holders, a couple of shoulder rests, handgrips, a paper roll holder and a stainless steel basin Ø 32 cm.

Delivered in kit form.

Made in Italy.