

# NEW OTO P/V ENT CHAIR - light blue



Product Code	27549
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2248674
CND	Z12149080
EAN/UPC	8023279275490

## Description

NEW OTO/PV ENT CHAIR with headrest - light blue Parigi

Patient chair with Peltex-lined seat and backrest.

Electric up and down movements by silent engine controlled through foot switch.

Back-rest can be adjusted up to  $-5^{\circ}$   $+90^{\circ}$  by a gas spring.

Chair rotates  $360^{\circ}$  and is provided with a rotation locking device placed at the bottom.

Power: 230 V - 50 Hz (Available 110 V - 60 Hz).

Max load 150 kg

Made in Italy.