

DONOR WHEELCHAIR - blue



|               |                |
|---------------|----------------|
| Product Code  | 27561          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |               |
|------------|---------------|
| UK-REP     | No            |
| CH-REP     | No            |
| RDM (NSIS) | 1229071       |
| CND        | Z129008       |
| EAN/UPC    | 8023279275612 |

Description

DONOR WHEELCHAIR

Donor wheelchair with structure made of epoxy painted oval and round steel tube.  
Backrest adjustable by gas pump.  
Padded leg rest manually adjustable.  
Metal footrest extractable. Armrest with articulations, manually adjustable in height.  
Padding made of high-density polyurethane foam, fireproof Class 1 IM.  
Covering made of washable, plastic material (leatherette), fireproof Class 1 IM.  
Back push handle.  
Mattress overall size: 166x55 cm.  
Carton size: 185x67xh 58 cm  
Four wheels Ø 100 mm, 2 directional and 2 with brake.  
Net weight: 53 kg.  
Gross weight: 55 kg.  
Maximum load: 150 kg.

Supplied with blanket tray with belts to fix oxygen bottle, IV stand and chart holder support.

Made in Italy.