

DONOR WHEELCHAIR with wheel - mechanical - white



|               |                |
|---------------|----------------|
| Product Code  | 27567          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |               |
|------------|---------------|
| UK-REP     | No            |
| CH-REP     | No            |
| RDM (NSIS) | 2255152       |
| CND        | V080202       |
| EAN/UPC    | 8055136920113 |

Description

DONOR ARMCHAIR WITH WHEELS - white

Strong and stylish armchair with structure made of oval, round and square metal tubing fully coated with white epoxy powder (legs oval tubing 50x25x1.5 cm; central frame square tubing 30x30x1.5 cm, cushions round tubing 20xØ 1.5 cm). High density polyurethane foam platform with wooden base covered by leatherette upholstery. Flame resistant materials Class 1 IM.  
Synchronized adjustment of the two cushions by a gas spring (cylinder/piston) allows a quick moving into TR/RTR position.  
Backrest inclination from +60° to -5° in Trendelenburg position, leg rest from -80° to +5°.  
Reclinable arm-rests can be height adjusted or rotated up to 180° when removed.  
4 castors, 2 with directional brake ensure mobility.  
Load: 150 kg.  
Size: 175x55xh 51 cm.

Multilingual manual: GB, FR, IT, ES, PT, DE, SE, PL, RO, GR

Made in Italy.