

STANDARD PLUS BED



| | |
|---------------|----------------|
| Product Code | 27607 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | No |
| CH-REP | Yes |

| | |
|------------|---------------|
| RDM (NSIS) | 2537730 |
| CND | V080602 |
| EAN/UPC | 8023279276077 |
| PARAF | 903113239 |
| GMDN | 34873 |

Description

STANDARD PLUS BED

Chrome-plated steel bed ends. Laminated plastic panels.
Rigid mesh mattress platform.
Adjustable back-rest.
Size: 206x89xh 90 cm.
Load: 110 Kg.

Made in Italy.