

STANDARD PLUS BED - with wheels 100 mm



| | |
|---------------|----------------|
| Product Code | 27608 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | No |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 2537730 |
| CND | V080602 |
| EAN/UPC | 8023279276084 |
| PARAF | 903113241 |
| GTIN | 08023279276084 |
| GMDN | 34873 |

Description

STANDARD PLUS BED with Ø 100 mm castors

Chrome-plated steel bed ends with Ø 100 mm castors. Laminated plastic panels.
Rigid mesh mattress platform.
Adjustable back-rest.
Size: 206x89xh 90 cm.
Load: 110 Kg.

Made in Italy.