



Product Code	27610
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2537730
CND	V080602
EAN/UPC	8023279276107
PARAF	903113254
GTIN	08023279276107
GMDN	34873

## Description

### BED WITH ADJUSTABLE HEAD AND FOOT SECTION

- BED WITH 3 ARTICULATIONS

Painted steel frame with net.  
 3 bed frames with nets in painted steel pipe, 200 x 85 cm.  
 2 Heads in chromed pipe complete with MDF bilaminated coated panel.  
 Movement by 2 cranks.  
 Protection of end run.  
 Rubber wheels of Ø 100 mm.  
 Size 206 x 89 x h 90 cm.  
 Height of frame from bottom 52/55 cm.  
 Size of frame: 200 x 85 cm.  
 Load: 110 Kg.  
 Weight 42 kg.

Made in Italy.