



Product Code	27610
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2537730
CND	V080602
EAN/UPC	8023279276107
PARAF	903113254
GMDN	34873

Description

BED WITH ADJUSTABLE HEAD AND FOOT SECTION

- BED WITH 3 ARTICULATIONS

Painted steel frame with net.
 3 bed frames with nets in painted steel pipe, 200 x 85 cm.
 2 Heads in chromed pipe complete with MDF bilaminated coated panel.
 Movement by 2 cranks.
 Protection of end run.
 Rubber wheels of Ø 100 mm.
 Size 206 x 89 x h 90 cm.
 Height of frame from bottom 52/55 cm.
 Size of frame: 200 x 85 cm.
 Load: 110 Kg.
 Weight 42 kg.

Made in Italy.