

KING HEIGHT ADJUSTABLE EXAMINATION COUCH - blue



Product Code	27613
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2262153
CND	V08060199
EAN/UPC	8023279276138

Description

KING HEIGHT ADJUSTABLE COUCH - blue - 9 - white frame

Height adjustable (48-95 cm) examination couch by electric actuator (Linak engine 6.000 N) controlled by foot switch with head section easily adjustable (to +60° and -10°) by a gas pump.

Strong and stylish leg structure in oval tubing 50x25x2 mm, central frame in rectangular tubing 50x20x1.5 mm painted with RAL 9003 white epoxy powder. Minimum thickness 70/80 micron, glossy finish.

High density polyurethane foam platform with wooden base covered by leatherette upholstery. Both materials are fireproof Class 1 IM.

Leatherette washable with the most common detergents.

Size: 190x60xh 48/95 cm.

Weight: 70 kg

Maximum load: 150 kg.

Operating power: 100-240 V, 50/60 Hz.

Couch roll holder included.

Made in Italy.

Technical Specifications

- Size: 190 x 60 x 48/95 cm
- Weight: 70 kg