



Product Code	27644
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2277577
CND	V0899
EAN/UPC	8023279276442
PARAF	903113355
GMDN	37510

## Description

### LIFT POLE

Lift pole made of steel frame with coating surface.  
The frame is composed of 4 parts, easy to assemble.  
Supplied with hand ring.

Multilingual manual: GB, FR, IT, ES, PT, DE, GR, RO, PL, SE.

## Technical Specifications

Overall length: 725 mm  
Overall width: 690 mm  
Overall height: 1935 mm  
Ring center: 1520 mm  
Load\*: 75 kg

\*Lifting stand load considers half body lift