



Product Code	27645
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2302554
CND	V0899
PARAF	987472646

Description

LIFTING STAND

Tubular painted steel frame with adjustable chromed steel joint.
Can be fully disassembled to reduce size.
Three main height adjustments.
Triangle hangs from belt with adjustable length.

Adjustable height 182 to 200 cm
Base 87x70 cm

Technical Specifications

Overall length: 864 mm
Overall width: 700 mm
Overall height: 2000 mm
Ring center: -
Load*: 80 kg
*Lifting stand load considers half body lift