

SPECIALIST BED 3 JOINTS/4 SECTIONS - electric with Trendelenburg



Product Code	27673
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	1368835
CND	Y181210
EAN/UPC	8023279276732
UMDNS	10-347
PARAF	987472370

Description

SPECIALISTIC BED with Trendelenburg - 40-80 cm

Three-joints, four sections bed with adjustable height from 40 to 80 cm with full electrical movements. All welded mainframe made of rectangular steel tube 50x25x2 mm. Mattress platform with flexible wooden staves, complete with mattress-hold guides. Backrest-section and leg-section adjustable by means of two independent electric motors 24V D.C. and push-button. Foot-section adjustable by means of a rack device. Fixed pelvis-section. Removable head and foot-end frames made of natural beech finish wood with solid corners and edges. Comes complete with 2 vertical sliding bars with upper/lower lock/unlock keys and lifting pole. Legs fitted with 4 castors Ø 100 mm with independent brake.

Technical Specifications

Size of platform: 900x2,000 mm
Platform height: 400-800 mm
External measurement: 1,030x2,140 mm
Max load: 135 kg
Weight: 102 kg
Operating power: 230 V - 50 Hz