

4 MOTORS, ELECTRICAL, HEIGHT ADJUSTABLE BED

24-65 cm - 3 joints - 4 sections



Product Code	27683
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2517670
CND	V08060101
EAN/UPC	8023279241716

Description

4 MOTORS ADJUSTABLE HEIGHT BED - 24-65 cm

Three-joints, four section beds with adjustable height by fully electrical movements.
Steel mainframe, surface with coating
Mattress platform with four sections, anti-sliding steel surface.
Backrest-section and leg-section adjustable by means of two independent electric motors.
Fixed pelvis-section.
Head and foot boards made of solid beech wood with solid corners and edges.
Legs fitted with 4 castors with independent brake.
Richmat 4 motors and 5 functions hand controller.
Supplied with lifting pole.

Multilingual manual: GB, FR, IT, ES, PT, SE, RO, GR.

Technical Specifications

Size of platform: 90x200 cm
Overall height: 74-114 cm
Platform height: 27-68 cm
External measurement: 102x215 cm
Max inclination of back section: 70°
Max inclination of leg section: 25°
Max inclination of overall bed: 12°
Height of bed end from mattress platform: 46 cm
Max load: 180 kg
Operating power: 230 V - 50 Hz