

TEKMA WATERPROOF COVER 170 g/m2 - white - universal



Product Code	27695
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2208409
CND	V0880
EAN/UPC	8023279276954
PARAF	987472077

Description

TEKMA WATERPROOF COVER 170 g/m2 - white - 195x85x18 cm

TEKMA fabric features excellent elasticity and softness, providing medical devices with high-performance qualities. It is flame-retardant, breathable, waterproof and antibacterial. Washable at 95°C and sterilizable up to 134°C.