



| | |
|---------------|----------------|
| Product Code | 27702 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| UK-REP | No |
| CH-REP | No |
| RDM (NSIS) | 2184713 |
| CND | Y122199 |
| EAN/UPC | 8023279277029 |
| PARAF | 903607137 |

Description

COMMODE WHEELCHAIR - painted

Powder coated steel frame, detachable armrest and footrest, solid rear castors with lock, plastic commode seat.

Technical Specifications

- Total width: 48 cm
- Seat width: 44 cm
- Rear wheel diameter: 12.5 cm
- Front wheel diameter: 12.5 cm
- Seat height: 51 cm
- Total wheelchair height: 92 cm
- Total length: 90 cm
- Seat depth: 44 cm
- Seat height: 42 cm
- Max. load: 100 kg
- Kit size: 54 x 53 x 53 cm
- NW GW (kg): 17-18