



Product Code	27704
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2692654
CND	Y1899
EAN/UPC	8023279277043

## Description

### FOLDING BEDSIDE RAIL

It can be attached to the bed to assist and help those people who have difficulty moving, getting up and lying down. One-click folding design frame: handle can be easily folded on both left and right side. Easy to use, single hand operation.

- stepless height adjustment
- 3 levels of handle width adjustment
- comfortable PU handle lowers and raises by 180 degrees
- equipped with light for night-time use
- made of iron

Size: 64-84x63.5xh 73.5-91.5 cm. Weight: 3.95 kg

Packaging size: 60x46x36 cm

Max load: 125 kg